

Remo Data Bridge User Manual

Version: 1.0 **Last Updated:** 2025-12-30

Table of Contents

1. [Accessing the Site](#)
 2. [Login](#)
 3. [Dashboard](#)
 4. [Request List](#)
 5. [Create New Request](#)
 6. [Request Detail](#)
-

1. Accessing the Site

Access URL

<http://remo-data-bridge.remo.re.kr:3015/>


How to Access

1. Open a web browser (Chrome, Edge, Safari, etc.).
2. Enter the URL above in the address bar.
3. Press Enter to access the site.

Recommended Browsers

- Google Chrome (latest version)
 - Microsoft Edge (latest version)
 - Safari (latest version)
-

2. Login



Description

- This is the login page for accessing the system.
- Enter your email and password to log in.


How to Use

1. Enter your email address.
2. Enter your password.
3. Click the **Login** button.

Account Creation

- If you do not have an account, you must request account creation from REMO.
 - Account creation request: info@remo.re.kr
-

3. Dashboard




Description

- This is the main screen displayed after logging in.
- You can view the overall request status at a glance.

Key Features

- View statistics by request status
 - View recent request list
 - Quick action menu
-

4. Request List



Description


- View all registered analysis requests in a list.

Key Features

- Search and filter requests
- Check request status (Pending, In Progress, Completed)
- Navigate to request detail page
- Create new request

5. Create New Request

5-1. Data Transfer Consent



Description

- When creating a new request, a **Cross-Border Data Transfer Consent** popup is displayed first.
- User consent is required as data will be transferred to overseas servers for analysis.

How to Use

1. Review the consent information.
 2. Click the consent checkbox at the bottom.
 3. Click the **Agree & Proceed** button.
 4. To decline, click the **Cancel** button to return to the previous page.
-

5-2. Enter Request Information

The screenshot shows the 'New Analysis Request' form in the REMO Data Bridge application. The form is titled 'New Analysis Request' and is located in the center of the screen. It has a dark blue header bar with the text 'REQUEST INFORMATION'. Below the header, there are three main sections: 'REQUEST TITLE' with a text input field containing 'ex) 2025-S-001', 'MEMO (OPTIONAL)' with a larger text area containing 'Enter additional notes or instructions', and 'UPLOAD FILES' with a dashed border box containing the text 'Drag & drop files or click to upload' and 'Maximum 10 files'. At the bottom of the form, there is a dark blue box with a shield icon and the text 'SECURITY DIRECTIVE' and 'AES-256 military-grade encryption active. Your data is protected during transfer and storage.' Below this box are two buttons: 'CANCEL' and 'CREATE REQUEST'.

Description

- This is the form for creating an analysis request after consent.


How to Use

1. Enter the **Request Title**.
2. **Upload files** required for analysis.
3. Add a **Memo** if you have additional notes.
4. Click the **Submit Request** button.

File Upload

- You can upload files by drag and drop.
 - Multiple files can be uploaded at once.
 - Uploaded files are securely transferred with AES-256 encryption.
-

6. Request Detail



Description

- This page displays detailed information about an individual analysis request.

Key Features

- View request information (title, status, creation date, etc.)
- View uploaded file list and download
- Check analysis results
- Download result files

File Download

- Click on a filename to download the file.

This document includes auto-generated screenshots.